

Title:	<b>SOP:</b> Covalent coupling of aSARS to NP
Employee:	Stefan Wagner, Surflay Nanotec GmbH
Protocol number:	FZS084, FZS080, FZS081

## 1 MATERIALS

- 0.1 % Particles 300 nm (Eu, QD, Blue, Green, Orange, Yellow)
  - deionized water (DI water)
  - PEG 20k
  - Tween20
  - anti-SARS 302 2 g/L (provided by FZMB)
  - EDC
  - S-NHS
  - MES buffer
  - Sodium acetate buffer (NaAc)
  - PBS
  - Ethanolamine
  - Sodium azide (azide)
- Required solutions:**
- Activation buffer (AB): 5 mM MES pH 5.5; 1 mM azide\*
  - Reaction buffer (RB): 5 mM NaAc pH 5.5; 0.5 % PEG 20k; 1 mM azide\*
  - Wash buffer (WB): 10 mM Hepes pH 7.4; 0.5 % PEG 20k; 0.01 % Tween20; 1 mM azide\*
  - 15 % PEG 20K
  - 0.3 % Tween 20
  - 10 mg/mL EDC in AB 5 mM MES (use immediately!)
  - 10 mg/mL S-NHS in AB 5 mM MES (use immediately!)
  - 25 % / 4 M Ethanolamine in DI water (use immediately!)

\*The Sodium azide is an optional component and is used for storing the buffer solutions in the fridge.

## 2 PARTICLE ACTIVATION

- 250  $\mu$ L particle suspension (0.1 %)
- centrifuge for 3-5 min at 12400xg \*\*; discard supernatant
- add 235  $\mu$ L AP, resuspend carefully and use Ultrasonic bath 10-30 s \*\*\*
- add 5  $\mu$ L EDC and 10  $\mu$ L S-NHS (c 10 mL/mL in AP)
- incubate 30 min at room temperature on a shaker 1600 rpm
- centrifuge for 3-5 min 12400xg \*\* and discard supernatant

## 3 COVALENT COUPLING

- resuspend the pellet carefully in 0,25 mL RB, use ultrasonic bath 100% 10-30 s, if necessary \*\*\*
- add 3.13-6.25  $\mu$ L of antibody (concentration 2 g/L, 3.13  $\mu$ L for research, 6.25  $\mu$ L for customers and QD)
- incubate for 60 min at room temperature on a shaker 1600 rpm
- add 5  $\mu$ L of 25 % Ethanolamine, 5  $\mu$ L 15 % PEG 20K and 5  $\mu$ L 0.3 % Tween20 to the mix
- incubate for an additional 30 min at room temperature on a shaker 1600 rpm
- centrifuge for 3-5 min at 12400g \*\*; discard supernatant
- wash twice with 250  $\mu$ L WB, use ultrasonic bath 100% 30 s, if necessary \*\*\*
- resuspend in 250  $\mu$ L WB (final concentration theor. 0.1 %)
- in the case of large losses due to washing, take up in less WB (e.g. 200  $\mu$ L) and determine the concentration by UV/Vis (10  $\mu$ L+190  $\mu$ L MilliQ in a micro cuvette)

\*\* Suggested centrifugation speed for Batch PSF007\_52, it can be further optimized to reduce agglomeration. Note: Centrifugation speed might change after activation and coupling.

\*\*\* Always check if it is at room temperature otherwise the Particles or the antibody can be damaged, ultrasonic bath is mostly optional but can improve performance